MARS 600 • 800 • 1000 • 1200 • 2000


Red light panel

User manual

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Mode d'emploi

Handleiding

Version 2.0 • 10/2026

RESTY-MARS-600-2000 • EN / DE / FR / NL

Choose your language

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Read the safety instructions before your first session and keep this manual with your panel.

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MARS 600 • 800 • 1000 • 1200 • 2000

CONTENTS

1	Welcome to your MARS panel	4
2	Safety	5
3	Your panel	6
4	Installation	7
5	Operating the touchscreen	8
6	Control with the app	10
7	Settings and effect	10
8	Care and storage	11
9	Troubleshooting	12
10	Technical data	13
11	Support, warranty and disposal	14

IMPORTANT

Read this manual carefully before using the panel and keep it for future reference. If you pass the panel on to someone else, include this manual.

SYMBOLS IN THIS MANUAL

WARNING

Risk of injury or damage. Always follow these instructions.

NOTE

Useful information for safe and correct use.

TIP

Practical advice.

1 Welcome to your MARS panel

Your Resty MARS panel combines nine LED wavelengths, from visible blue, yellow and red light to invisible near-infrared. This manual applies to the MARS 600, 800, 1000, 1200 and 2000. All five models share the same touchscreen, app and settings. They differ only in size, number of LEDs and power (see Technical data).

Quick start

- 1 Check the panel, power cable, protective glasses and remote control for damage.
- 2 Hang the panel on a door with the supplied hanging kit (see Installation), or mount it on a Resty stand according to the instructions supplied with the stand.
- 3 Plug the power cable into the socket on the back of the panel and into a wall outlet. Switch the panel on with the power switch on the back.
- 4 Put on the protective glasses.
- 5 Tap **Scene** to choose a preset, or set the timer and mode under **Settings**. Tap **Start/Pause** to start.
- 6 Keep a distance of 15 to 30 cm. Stop immediately if anything feels uncomfortable.

NOTE

General wellness use. This manual does not provide treatment protocols for medical conditions and does not replace medical advice. Never use the panel while sleeping.

NOTE

Screen labels. The panel's screen is available in English, German and Chinese. This manual shows the screen labels in English, for example **Scene**, **Settings** and **Start/Pause**. The labels on your screen may differ slightly.

2 Safety

Read all warnings before you use the panel. Ignoring them can result in injury, electric shock or damage to the device.

Light and eyes

- Wear the supplied protective glasses during every session. Do not look directly into the LEDs while they are on.
- Near-infrared light (NIR) is largely invisible. LEDs that look dark can still emit light.

Do not use the panel if you

- are sensitive to light or take medication that makes you more sensitive to light;
- have ever had a seizure triggered by light;
- have skin type V or VI on the Fitzpatrick scale;
- are pregnant, breastfeeding or planning a pregnancy;
- have had sunburn in the past two weeks;
- have open wounds in the area you want to expose.

During and after use

- Stop immediately if you feel any discomfort. Consult a doctor if necessary.
- If you suspect an adverse reaction, stop using the panel and seek medical advice.
- Children and persons who cannot recognise discomfort or end a session on their own must not use the panel independently. Never use the panel on a person who is asleep or unconscious.

Electrical safety

This is an electrical device. To prevent electric shock and fire:

- Never let the panel come into contact with water or other liquids, and never clean electrical parts with liquids.
- Do not touch the panel, plug or switches with wet hands.
- Do not pull, carry or lift the panel by its cable.
- Do not use the panel if the cable is damaged, or if the panel has been dropped, left outdoors or exposed to water or other liquids.
- Unplug by holding the plug, never by pulling the cable.
- Unplug the panel when it is not in use, and route the cable so that nobody can trip over it.
- Do not use the panel near flammable or combustible products, or where vapours may be present: risk of fire or explosion.
- Do not use the panel on or near heated surfaces.

Device safety

- If the panel overheats, switch it off immediately.
- Keep the ventilation openings and cooling elements on the back clear. Do not cover the panel during use.
- Do not open, disassemble or repair the panel yourself.
- Only use accessories supplied or specified by the manufacturer. Other accessories can increase electromagnetic emissions or reduce electromagnetic immunity.
- Do not scratch, damage or drop the panel. For indoor use only.
- Store the panel in a cool, dry place away from direct sunlight and hot surfaces. High temperatures can damage the device.
- Keep cables out of the reach of children: risk of strangulation.
- Always operate, transport and store the panel as described in this manual. Failure to do so voids the warranty.

WARNING

Burning smell, smoke, unusual heat or a damaged cable? Switch the panel off, unplug it if this is safe, and contact Resty before you use it again.

3 Your panel

FRONT

BACK

Illustrations show an example; your model may differ in size.

LED surface

Front of the panel: nine wavelengths behind 30° lenses.

Touchscreen	Scenes, timer and all settings.
Power switch	On the back. Switches the power supply on and off.
Power socket	On the back. For the supplied power cable.
Cascade connection	On the back, for linking panels. Only use it with original Resty cables and after confirmation from Resty.
Cooling elements	On the back. Keep them clear during use.

In the box

Panel, power cable, protective glasses, remote control, this manual and, depending on the model, a door-hanging kit with threaded posts, braided hanging cables, a door hook and a height-adjustable pulley. Contact Resty if anything is missing or damaged.

Remote control

Press ►|| on the remote control to start or pause the light. With **SET**, **+** and **–** you can change settings from a distance; the small display on the remote control shows what you are adjusting. If the remote control does not respond, check its battery.

4 Installation

Hanging the panel on a door

- 1 Place the panel face down on a clean, scratch-free surface. Screw the threaded posts into the top of the panel and attach the braided cables securely to the posts.
- 2 Hang the door hook over the top of a sturdy, standard-sized door.
- 3 Connect the braided cables to the snap hook to hang the panel on the door.
- 4 To raise the panel, support it with one hand and pull down the loose end of the cable with the other hand.
- 5 To lower the panel, push up the lever of the height-adjustable pulley while supporting the panel with your other hand. Always keep one hand on the panel to support its weight while adjusting.
- 6 Plug the power cable into the socket on the back of the panel and the other end into a wall outlet. Switch the panel on with the power switch.

TIP

Use a door that stays closed during the session, and ask a second person to help with the larger models.

On a stand

Resty stands come with their own assembly instructions. Always follow those instructions and respect the maximum load of the stand.

5 Operating the touchscreen

TOUCHSCREEN

Main screen

Display	Shows the timer and mode you have set.
Scene	Choose a preset: Fat Loss, Skin, Muscle, Joint Pain, Sleep, Pets, Universal or Custom. See Settings and effect.
Start/Pause	Switches the light on or pauses it.
Settings	Timer, mode, brightness, frequency, group, alarm clock, date and time, and other settings.

1 • Timer

Adjustable from 1 to 30 minutes. Drag the marker along the circle, or tap the centre of the circle, enter the time and confirm with the check mark. If you do not set a time, the panel switches off after 5 minutes.

2 • Mode

Choose which wavelengths are on. Wavelengths in nanometres (nm):

ALL	480 590 630 660 670 810 830 850 1060
RED&NIR&Yellow	590 630 660 670 810 830 850 1060
RED&Yellow	590 630 660 670
RED&NIR	630 660 670 810 830 850 1060
NIR&Yellow	590 810 830 850 1060

RED	 630  660  670
Yellow	 590
NIR	 810  830  850  1060
Blue	 480

3 • Brightness

Adjust the brightness of each wavelength separately from 0 to 100 %: tap **–** or **+**, use the slider, or tap the % value, enter a number and confirm with the check mark.

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Set a pulse frequency for each wavelength from 0 to 10,000 Hz: tap **–** or **+**, use the slider, or tap the value and enter a frequency. If no frequency is set, the light shines continuously.

5 • Group

Set one panel as **Primary** and the other panels as **Secondary**. All panels are then controlled by the primary panel. Tap **Close** to end group control.

6 • Alarm clock

- Choose **Only once** or **Repeat every day**, then set the alarm time, timer, mode, brightness and pulse.
- The alarm time is the moment the session ends. Example: with a 10-minute timer, the light switches on 10 minutes before the alarm time.
- Swipe left to switch an alarm off, swipe right to switch it on.
- Two alarms are available by default. Add or delete alarms with **+** or **–**, up to 10.

7 • Date & time

Set the date and time for your region.

8 • Other settings

- **Bluetooth:** switch on to connect the app.
- **Work alarm:** the panel beeps when half of the session time has passed, as a reminder to turn to another body area.
- **Volume** of the operating sounds: 0, 25, 50, 75 or 100 %.
- **Screen brightness:** 0, 25, 50, 75 or 100 %.
- **Factory reset:** restores the default settings.
- **Language:** English, Deutsch, 简体中文 or 繁體中文.

9 • About

Shows the device information, including the Bluetooth name (BLE name) you need for the app.

6 Control with the app

Instead of the touchscreen you can use the free app. It offers the same functions: timer, mode, brightness, pulse and start/pause.

iPhone / iPad App Store: search for **IDEA-BLE**

Android Google Play: search for **IDEA** (provider: Red light therapy)

- 1 On the panel, switch on Bluetooth: **Settings → Other Setting → Bluetooth**.
- 2 Switch on Bluetooth on your phone. Allow the app to use Bluetooth when it asks.
- 3 Open the app and connect to your panel. You find the panel's Bluetooth name (BLE name) under **Settings → About**.

NOTE

The app is in English. Keep your phone within a few metres of the panel while connecting.

7 Settings and effect

What the wavelengths do

Visible light 480 nm blue • 590 nm yellow (amber) • 630, 660 and 670 nm red. Often chosen for complexion care and skin-focused routines.

Near-infrared 810, 830, 850 and 1060 nm. Invisible to the eye and reaches deeper than visible red light. Often included in muscle and recovery-focused routines.

Which mode for which routine?

Routine	Suggested mode
Skin and face	RED or RED&Yellow; Blue can be added for skin-focused sessions
After training or a long day	RED&NIR or NIR
All-round routine	ALL (all nine wavelengths) or RED&NIR

Evening wind-down

RED; avoid Blue in the evening

Distance, time and frequency

- **Distance:** 15 to 30 cm from the LED surface to your skin. Start at about 30 cm.
- **Session:** start with 10 minutes per body area; 10 to 15 minutes is a common session length. The timer goes up to 30 minutes.
- **How often:** a few times per week, for example three times to start. Regular use matters more than long sessions.
- **Rest:** wait at least 6 hours between two 30-minute sessions.
- **Brightness:** if the light feels too intense or warm, increase the distance or lower the brightness. Do not compensate with a longer session.
- **Pulse:** optional. If you are unsure, use continuous light (no frequency set).
- Use the panel on clean, dry, uncovered skin, and keep the same distance throughout the session.

Scenes

The scene names (Fat Loss, Skin, Muscle, Joint Pain, Sleep, Pets, Universal and Custom) come from the panel software. Each scene loads a preset combination of settings. The main screen shows the active timer and mode, and you can change them under **Settings**. The scenes are not medical treatments.

NOTE

MARS panels are designed for general wellness, skin-care and fitness routines. They are not intended to diagnose, treat, cure or prevent any disease.

8 Care and storage

- Switch the panel off, unplug it and let it cool down before cleaning.
- Wipe the housing and LED surface with a soft, dry cloth. Do not use sprays, solvents or abrasive cleaners.
- Check the panel, cable and plug before each use. Do not use damaged parts.
- Store the panel in a cool, dry place, away from direct sunlight and heat.
- Replace the protective glasses if they are scratched or damaged.

9 Troubleshooting

Problem	Possible cause	Solution
The panel does not work	Not plugged in	Plug the power cable into the socket on the back and into a wall outlet.
	Not switched on	Press the red power switch on the back.
The panel suddenly stops	The set session time has ended. Without a set time, the panel switches off after 5 minutes.	Start a new session if you wish. Wait at least 6 hours between two 30-minute sessions.
No pulsed light	No frequency set	Set a frequency (0 to 10,000 Hz) under Settings with + or –.
The LEDs seem to be off	Near-infrared is invisible to the eye	Check the mode. Light is visible in modes with red, yellow or blue. Do not look into the LEDs to check.
The app does not find the panel	Bluetooth is off or app permission is missing	Switch on Bluetooth on the panel and on your phone, allow Bluetooth access and check the BLE name under Settings → About.
The remote control does not respond	Battery empty	Check or replace the battery of the remote control.

WARNING

Do not attempt to repair the panel yourself. Contact Resty if a problem persists.

10 Technical data

	MARS 600	MARS 800	MARS 1000	MARS 1200	MARS 2000
Dimensions (cm)	32 × 23 × 6.5	64 × 23 × 6.5	92 × 23 × 6.5	92 × 42 × 7	180 × 41 × 6.5
LEDs	72	144	216	432	864
Weight (kg)	3.9 ± 0.5	6.0 ± 0.5	8.5 ± 0.5	15.3 ± 0.5	29.4 ± 0.5
Rated power (W)	120	240	360	720	1440
Irradiance at 15 cm, R+NIR (mW/cm²)	157.6	163.6	175.1	181.3	≥ 140
Current at 220 V (A)	0.918	1.735	2.561	2.994	12
Current at 110 V (A)	1.549	3.019	4.401	6.113	16
Wavelengths	480, 590, 630, 660, 670, 810, 830, 850 and 1060 nm; 9 modes				
Input voltage	AC 100—240 V				
Controls	Touchscreen, remote control and app (Bluetooth)				
Brightness / pulse	Per wavelength 0—100 % / 0—10,000 Hz				
Timer	1—30 minutes				
Lens angle	30°				
EMF	0.0 µT at 10 cm				
Temperatures	Operation 5—50 °C • storage -20 to +60 °C • surface ≤ 40 °C				
Rated LED lifespan	50,000 hours				
Housing	Metal				

Current values ± 5 %. Irradiance and current are manufacturer specifications.

11 Support, warranty and disposal

Contact	info@feelresty.com • feelresty.com
----------------	------------------------------------

Warranty	Bought on feelresty.com: 3-year warranty from the date of receipt. Bought elsewhere: the warranty terms of your seller apply. Your statutory rights are not affected. Opening or modifying the panel voids the warranty.
-----------------	--

Disposal	Do not dispose of the panel with household waste. Take it to a collection point for electrical and electronic equipment.
-----------------	--

When you contact us, please include your order number, the model and a short description of the problem. A photo or video often helps.

Reventis B.V. • brand Resty • Vlierweg 12, K 03 • 1032 LG Amsterdam • The Netherlands

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Reventis B.V. • merk Resty • Vlierweg 12, K 03 • 1032 LG Amsterdam • Nederland
Reventis B.V. • Resty